

Dilated Cardiomyopathy/Left Ventricular Noncompaction Panel

Disorder also known as: Idiopathic Dilated Cardiomyopathy (IDC); Familial Dilated Cardiomyopathy (FDC); Left Ventricular Noncompaction Cardiomyopathy (LVNC)

Panel Gene List: *ABCC9, ACTC1, ACTN2, ALMS1, ANKRD1, BAG3, CHRM2, CRYAB, CSRP3, DES, DMD, DOLK, DSC2, DSG2, DSP, DTNA, EMD, FKTN, GATAD1, ILK, LAMA4, LAMP2, LDB3, LMNA, MIB1, MTND1, MTND5, MTND6, MTTD, MTTT, MTTI, MTTK, MTTL1, MTTL2, MTTM, MTTQ, MTTTS1, MTTTS2, MYH6, MYH7, MYBPC3, MYPN, NEBL, NEXN, PLN, PRDM16, RAF1, RBM20, SCN5A, SGCD, TAZ, TCAP, TMPO, TNNC1, TNNI3, TNNT2, TPM1, TTN, TTR, TXNRD2, VCL*

Additional genes from our cardiology test menu may be added to this panel by selecting test code J554C.

Clinical Features:

Dilated cardiomyopathy (DCM) is a disease of the heart muscle that is diagnosed based on the findings of both left ventricular enlargement and systolic dysfunction, resulting in a reduction in the myocardial force of contraction.^{1,2} This condition is primarily determined by echocardiogram to measure cardiac chamber dimensions, ventricular thickness, and ejection fraction.^{2,3} DCM usually presents with one or more of the following; i) heart failure with symptoms of congestion (edema, orthopnea or paroxysmal dyspnea), ii) reduced cardiac output, resulting in fatigue or dyspnea on exertion, arrhythmias and/or conduction system disease and iii) thromboembolic disease or stroke, mainly from left ventricular mural thrombus.² However, individuals with DCM may also be asymptomatic. The prevalence of idiopathic DCM is at least 1/2,700 in the general population.²

Left ventricular non-compaction (LVNC) is characterized by abnormal trabeculations in the left ventricle, most frequently at the apex. Prominent trabeculations and deep intertrabecular recesses give the myocardium a sponge-like appearance that can be detected using echocardiogram.^{1,3,4} LVNC can share the same clinical presentation as DCM, ranging from asymptomatic disease to progressive deterioration of cardiac function, arrhythmias, thromboembolic events, or sudden cardiac death.^{3,4,5} LVNC is a rare condition, which affects less than 0.3% of the population.³

While DCM/LVNC can be an isolated finding, it can also occur in association with neuromuscular disorders or mitochondrial myopathy.^{2,4,5} In addition, DCM/LVNC can be a presenting feature of various genetic syndromes, including Danon disease, Carvajal Syndrome, Barth syndrome, or Emery-Dreifuss muscular dystrophy.^{2,4,5}

Inheritance Pattern/Genetics: Autosomal Dominant, Autosomal Recessive, X-linked, or Mitochondrial

Test Methods:

Using genomic DNA from the submitted specimen, the coding regions and splice junctions of the 61 genes (except the following genomic regions of the *TTN* gene: chr2:179527692-179527782, 179523898-179523982, 179523731-179523815) are enriched using a proprietary targeted capture

system developed by GeneDx. These targeted regions are sequenced simultaneously by massively parallel (NextGen) sequencing on an Illumina platform with paired-end reads. Bi-directional sequence is assembled, aligned to reference gene sequences based on human genome build GRCh37/UCSC hg19, and analyzed for sequence variants. Capillary sequencing is used to confirm all potentially pathogenic variants and to obtain sequence for regions where fewer than 15 reads are achieved by NextGen sequencing. Concurrent deletion/duplication testing is performed for the genes in the panel using exon-level oligo array CGH (ExonArrayDx), except for the 13 mitochondrial genes. The *EMD* and *TAZ* genes have gene level resolution; exon level events may not be detected. Data analysis is performed using gene-specific filtering. Probe sequences and locations are based on human genome build GRCh37/UCSC hg19. The array is designed to detect most intragenic deletions and duplications. Confirmation of copy number changes is performed by MLPA, qPCR, or repeat array CGH analysis. Sequence and array CGH alterations are reported according to the Human Genome Variation Society (HGVS) or International System for Human Cytogenetic Nomenclature (ISCN) guidelines, respectively. Benign and likely benign variants, if present, are not included in this report but are available upon request.

Test Sensitivity:

The clinical sensitivity of sequencing and deletion/duplication analysis of the 61 genes included in the DCM/LVNC Panel depends in part on the patient's clinical phenotype and family history. In general, the sensitivity is highest for individuals with clearly defined DCM/LVNC and a family history of disease. The technical sensitivity of the sequencing test is estimated to be 98%. The sequencing panel will not reliably detect deletions, insertions, or rearrangements greater than or equal to five base pairs (bp). Deletions or duplications of less than 500 bp are not reliably detected by array CGH.

Gene	Protein	Inheritance	Disease Association(s)
<i>ABCC9</i>	ATP-BINDING CASSETTE, SUBFAMILY C, MEMBER 9	AD	DCM, Cantu syndrome
<i>ACTC1</i>	ACTIN, ALPHA, CARDIAC MUSCLE	AD	CHD, DCM, HCM, LVNC
<i>ACTN2</i>	ACTININ, ALPHA-2	AD	Pediatric HCM/DCM
<i>ALMS1</i>	CENTROSOME AND BASAL BODY ASSOCIATED PROTEIN	AR	Alstrom syndrome, mitogenic cardiomyopathy
<i>ANKRD1</i>	ANKYRIN REPEAT DOMAIN-CONTAINING PROTEIN 1	AD	HCM, DCM
<i>CHRM2</i>	M2-MUSCARINIC ACETYLCHOLINE RECEPTOR	AD	DCM
<i>CRYAB</i>	CRYSTALLIN, ALPHA-B	AD, AR	DCM, myofibrillar myopathy
<i>CSRP3</i>	CYSTEINE- AND GLYCINE-RICH PROTEIN 3	AD	HCM, DCM
<i>DES</i>	DESMIN	AD	DCM, ARVC, myopathy, AV block
<i>DMD</i>	DYSTROPHIN	XL	DMD, BMD, DCM
<i>DOLK</i>	DOLICHOL KINASE	AR	DCM, congenital disorder of glycosylation type Im
<i>DSC2</i>	DESMOCOLLIN	AD, AR	ARVC, ARVC+skin and hair findings
<i>DSG2</i>	DESMOGLEIN	AD	ARVC, DCM
<i>DSP</i>	DESMOPLAKIN	AD, AR	ARVC, DCM, Carvajal syndrome
<i>DTNA</i>	DYSTROBREVIN, ALPHA	AD	LVNC, CHD
<i>EMD</i>	EMERIN	XL	EMD
<i>FKTN</i>	FUKUTIN	AR	DCM, LGMD, Fukuyama Congenital Muscular Dystrophy
<i>GATAD1</i>	GATA ZINC FINGER DOMAIN-CONTAINING PROTEIN 1	AR	DCM
<i>ILK</i>	INTEGRIN-LINKED KINASE	AD	DCM
<i>LAMA4</i>	LAMININ, ALPHA-4	AD	DCM

<i>LAMP2</i>	LYSOSOME-ASSOCIATED MEMBRANE PROTEIN 2	XL	Danon disease
<i>LDB3</i>	LIM DOMAIN-BINDING 3	AD	DCM, LVNC, myopathy
Gene (cont.)	Protein	Inheritance	Disease Association(s)
<i>LMNA</i>	LAMIN A/C	AD, AR	DCM, HCM, congenital muscular dystrophy, EMD
<i>MIB1</i>	MINDBOMB E3 UBIQUITIN PROTEIN LIGASE 1	AD	LVNC
<i>MTND1</i>	mtDNA ENCODED COMPLEX I, SUBUNIT ND1	MITO	Cardiomyopathy, myopathy
<i>MTND5</i>	mtDNA ENCODED COMPLEX I, SUBUNIT ND5	MITO	Cardiomyopathy, myopathy
<i>MTND6</i>	mtDNA ENCODED COMPLEX I, SUBUNIT ND6	MITO	Cardiomyopathy, myopathy
<i>MTTD</i>	MITOCHONDRIAL tRNA FOR ASPARTIC ACID	MITO	Cardiomyopathy, myopathy
<i>MTTH</i>	MITOCHONDRIAL tRNA FOR HISTIDINE	MITO	Cardiomyopathy, myopathy
<i>MTTI</i>	MITOCHONDRIAL tRNA FOR ISOLEUCINE	MITO	Cardiomyopathy, myopathy
<i>MTTK</i>	MITOCHONDRIAL tRNA FOR LYSINE	MITO	Cardiomyopathy, myopathy
<i>MTTL1</i>	MITOCHONDRIAL tRNA FOR LEUCINE 1	MITO	Cardiomyopathy, myopathy
<i>MTTL2</i>	MITOCHONDRIAL tRNA FOR LEUCINE 2	MITO	Cardiomyopathy, myopathy
<i>MTTM</i>	MITOCHONDRIAL tRNA FOR METHIONINE	MITO	Cardiomyopathy, myopathy
<i>MTTQ</i>	MITOCHONDRIAL tRNA FOR GLUTAMINE	MITO	Cardiomyopathy, myopathy
<i>MTTS1</i>	MITOCHONDRIAL tRNA FOR SERINE 1	MITO	Cardiomyopathy, myopathy
<i>MTTS2</i>	MITOCHONDRIAL tRNA FOR SERINE 2	MITO	Cardiomyopathy, myopathy
<i>MYBPC3</i>	MYOSIN-BINDING PROTEIN C, CARDIAC	AD	HCM, DCM
<i>MYH6</i>	MYOSIN, HEAVY CHAIN 6, CARDIAC MUSCLE, ALPHA	AD	CHD, DCM, HCM, SSS
<i>MYH7</i>	MYOSIN, HEAVY CHAIN 7, CARDIAC MUSCLE, BETA	AD	DCM, HCM, myopathy
<i>MYPN</i>	MYOPALLADIN	AD	DCM, RCM, HCM
<i>NEBL</i>	NEBULETTE	AD	DCM, endocardial fibroelastosis
<i>NEXN</i>	NEXILIN	AD	DCM, HCM
<i>PLN</i>	PHOSPHOLAMBAN	AD	DCM, HCM
<i>PRDM16</i>	PR DOMAIN CONTAINING 16	AD	DCM, LVNC
<i>RAF1</i>	V-RAF-1 MURINE LEUKEMIA VIRAL ONCOGENE HOMOLOG 1	AD	Noonan/CFC/Costello syndromes
<i>RBM20</i>	RNA-BINDING MOTIF PROTEIN 20	AD	DCM
<i>SCN5A</i>	SODIUM CHANNEL, VOLTAGE-GATED, TYPE V, ALPHA SUBUNIT	AD	BrS, DCM, Heart block, LQTS, SSS, SIDS
<i>SGCD</i>	SARCOGLYCAN, DELTA	AD, AR	DCM, LGMD
<i>TAZ</i>	TAFAZZIN	XL	DCM, LVNC, Barth syndrome
<i>TCAP</i>	TITIN-CAP (TELETHONIN)	AD, AR	HCM, DCM, LGMD
<i>TMPO</i>	THYMOPOIETIN	AD	DCM
<i>TNNC1</i>	TROPONIN C, SLOW	AD	DCM, HCM
<i>TNNI3</i>	TROPONIN I, CARDIAC	AD, AR	DCM, HCM, RCM
<i>TNNT2</i>	TROPONIN T2, CARDIAC	AD	DCM, HCM, RCM, LVNC
<i>TPM1</i>	TROPOMYOSIN 1	AD	DCM, HCM
<i>TTN</i>	TITIN	AD	DCM, myopathy
<i>TTR</i>	TRANSTHYRETIN	AD	TTR-related amyloidosis
<i>TXNRD2</i>	THIOREDOXIN REDUCTASE 2	AD, AR	DCM
<i>VCL</i>	VINCULIN	AD	HCM, DCM, LVNC

Abbreviations: AD – Autosomal dominant; AF – Atrial fibrillation; AR – Autosomal recessive; ARVC – Arrhythmogenic Right Ventricular Cardiomyopathy; AV block- Atrioventricular Block; BMD – Becker Muscular Dystrophy; BrS – Brugada Syndrome; CCD- Cardiac Conduction Disease; CHD – Congenital Heart Defects; CPVT – Catecholaminergic Polymorphic Ventricular Tachycardia; DCM – Dilated Cardiomyopathy; Abbreviations (cont.): DMD- Duchenne Muscular Dystrophy; EMD – Emery Dreifuss Muscular Dystrophy; HCM – Hypertrophic Cardiomyopathy; JLNS – Jervell and Lange-Nielsen Syndrome; LGMD – Limb Girdle Muscular Dystrophy; LQTS – Long QT Syndrome Abbreviations (cont.): LVNC – Left Ventricular Non-Compaction; RCM – Restrictive Cardiomyopathy; SIDS – Sudden Infant Death Syndrome; SSS – Sick Sinus Syndrome; TAAO- Thoracic Aortic Aneurysm and Dissection; XL – X-linked

References:

1. Online Mendelian Inheritance in Man, OMIM. Johns Hopkins University, Baltimore, MD. MIM Number: {115200; 604169}: {04/21/2015; 09/04/2013}: World Wide Web URL: <http://omim.org/>
2. Hershberger RE, Morales A. Dilated Cardiomyopathy Overview. 2007 Jul 27 [Updated 2013 May 9]. In: Pagon RA, Adam MP, Ardinger HH, et al., editors. GeneReviews [Internet]. Seattle (WA): University of Washington, Seattle; 1993-2015. Available from: <http://www.ncbi.nlm.nih.gov/books/NBK1309/>
3. Callis et al. (2010) Expert review of molecular diagnostics 10(3):329-351 (PMID: 20370590)

4. Bennett et al. (2016) *Cardiology Research and Practice* 2016:5172308 (PMID: 26881173)
5. Rooms et al. (2015) *Acta Cardiol* 70 (6):625-31 (PMID: 26717209)